


## hCG Screening and Pregnancy Testing

# CLINITEST® hCG Pregnancy Test: clinical evidence

siemens-healthineers.com/poc-specialty-urine-testing

## Assessment of clinical performance and sensitivity

The acceptable clinical performance and sensitivity of the CLINITEST hCG test reaffirm its utility as a **rapid and reliable** method for detection of pregnancy at an early stage.<sup>1</sup>



- 100% of tests on all positive-spiked samples reported positive.
- 100% of tests on all negative samples reported negative.



Concordance with all sensitivity acceptance criteria

The incidence of false-positive results was 3 in 13,500 negative urine samples tested ( $\leq 0.04\%$ ).



An acceptably low incidence of false positives with urine samples negative for hCG

Results from 4725 individual tests showed a rate of agreement for positive results well above the acceptable criteria of  $\geq 99.1\%$ .



Excellent rate of agreement of positive results with hCG-positive samples

### Potential causes of false-positive results:<sup>1</sup>

- Undetected early spontaneous abortion
- The presence of interfering substances such as endogenous proteins, drugs, microorganisms, and erythrocytes/leukocytes
- The detection of low concentrations of hyperglycosylated hCG and beta core fragments

## Greater efficiency. Less error. Improved documentation.



When interfaced with the EMR, the CLINITEK Status® Connect System would eliminate most pre- and post-analytical errors, especially transcription errors and delayed/improper/lack of documentation of some test results.<sup>2</sup>



**Test time + total time**  
Semiautomated testing < manual



**Transcription errors**  
0.3–1.7% manual vs. 0% semiautomated



**Lack of documentation**  
20.5% lower than manual entry


## Assessment of false-negative risk due to high-dose hook interference sensitivity

The CLINITEST hCG Test-CLINITEK Status hCG device is **unlikely to exhibit false-negative urinary hCG results due to high-dose hook interference for women in early healthy pregnancy.**<sup>3</sup>


The **hook effect** occurs when the hCG concentration is so high that it saturates both antibodies and there are so many molecules that the antibodies do not actually form a sandwich, producing a false-negative result



### Max. Lattice Formation and Max. Precipitates



## AACC recommendations to reduce the hCG hook effect



- When clinically appropriate, **quantitative hCG measurement should be performed** in serum or plasma specimens to avoid hCG $\beta$ cf interference.
- If urine testing is required, a device with minimal susceptibility to hCG $\beta$ cf variant effect should be used.

If a false-negative result due to the hook effect is suspected:



- The urine specimen **be diluted** using an hCG-free material and retested.
- A **positive result** upon dilution confirms the hook effect.
- Alternatively, a serum hCG test can be conducted.<sup>4</sup>

### References:

1. CLINITEST hCG testing on the CLINITEK Status analyzer achieves clinical performance and sensitivity criteria: an internal validation study. Siemens Healthcare Diagnostics Inc. 2014. Order No. A91DX-POC-131828-GC1-4A00.
2. Young PE, Diaz GJ, Kalariya RN, Mann PA, Benbrook MC, Avandsalehi KR, Petersen JR. Comparison of the time required for manual (visually read) and semiautomated POCT urinalysis and pregnancy testing with associated electronic medical record (EMR) transcription errors. Clin Chim Acta. 2020;504:60-3.
3. Milhorn D, Korpi-Steiner N. Using a simulation model to assess risk of false negative point-of-care urinary human chorionic gonadotropin device results due to high-dose hook interference. Clin Biochem. 2015 Feb;48(3):99-104.
4. Nichols JH, Ali M, Anetor JJ, Chen LS, Chen Y, Collins S, et al. AACC guidance document on the use of point-of-care testing in fertility and reproduction. J Appl Lab Med. 2022;7:1202-36.

CLINITEK Status, CLINITEST, and all associated marks are trademarks owned by or licensed to Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners. Product availability may vary from country to country and is subject to varying requirements.